



300 Breesport, San Antonio, TX 78216

**FIRST, READ ALL THE WARNINGS
AND INSTRUCTIONS ON TCP CAN**

Then calculate the amount of TCP
required for each fuel tank. 30ML treats
10 US Gallons Blue 100LL Avgas.

Using the Syringe Kit – P/N: 73138

Remove the cap from the can. Pour out
the desired amount of TCP into the
clean polypropylene beaker.

Using syringe, pull back on the plunger
and draw the TCP from the beaker.

Carefully dispense/ inject TCP into
aircraft fuel tank taking care not to get
any on painted surfaces.

Repeat the process until all the desired
amount of TCP has been added to the
fuel tank.

To Use ALCOR TCP Fuel Treatment

1-800-354-7233

www.alcoraero.com